

Automated Fully Enclosed qPCR Instrument

Galaxy Neo

Core Technology

1

7 colors Fluorescence channel

3

User Friendly Software



Galaxy Nano

2

Extraction, Purification,
Amplification and Results Analysis
integrated in one system

4


Processes 2 samples simultaneously
Melting Curve Analysis
LIS/HIS support

Ready to use

Pre-loaded reagent
Automated lyophilization technology

Fully Enclosed

Contamination free, Biosafety guaranteed
No PCR lab required



Open System

Built-in Sample Pre-treatment
with Ultrasound Heating and Lysis

Compatibility


Extract nucleic acids from
complex biological samples

iCassette

Workflow



Add sample to the iCassette



load iCassette to the instrument



Result report automatically

Specifications

Heat ramping		≥15°C/s						
Fluorescent linear		R≥0.99						
Precision		CV≤3%						
Size		448*178*331mm						
Weight		9kg						
Fluorescence channel	Channel No.	1	2	3	4	5	6	7
	Excitation Filter	460nm	525nm	543nm	571nm	624nm	675nm	350nm
	Emission Filter	525nm	564nm	584nm	624nm	675nm	710nm	460nm
	Suitable Dyes	FAM	HEX	TAMRA	ROX/TEXAS RED	CY5	CY5.5	AMCA
	Light Source	white LED						UV LED

Ordering

Classification	Cat.No.	Product Name	Packaging	Multiplex Detection
Instrument	101-0062-02	Automated Fully Enclosed qPCR Instrument (Galaxy Neo)	1 unit	2 sample throughputs
Universal kit	106-0051-01	Nucleic Acid Extraction and Purification Kit (One in one)	12 tests/kit	up to 18-plexes
Respiratory Infectious Diseases	106-0057-01	SARS-2-CoV (Covid-19)/ Influenza A Virus/ Influenza B Virus/ Respiratory Syncytial Virus Diagnostic Kit	12 tests/kit	5-plexes
	106-0100-01	13 multiplex Respiratory Pathogens Detection Kit	12 tests/kit	13-plexes
	106-0052-01	Influenza A Virus/ Influenza B Virus/ Respiratory Syncytial Virus Diagnostic Kit	12 tests/kit	3-plexes
	106-0056-01	Diagnostic Kit for Mycobacterium Tuberculosis(MTB) complex DNA and Rifampicin resistance Mycobacterium tuberculosis	12 tests/kit	6-plexes
	106-0069-01	Mycoplasma Pneumoniae/ Adenovirus/ Rhinovirus/ Human Metapneumovirus Diagnostic Kit	12 tests/kit	4-plexes
	106-0055-01	Streptococcus Pneumoniae/ Klebsiella Pneumoniae/ Legionella Pneumophila/ Haemophilus influenzae/ Bordetella Pertussis Diagnostic Kit	12 tests/kit	5-plexes
Hospital Acquired Infection (HAI)	106-0102-01	Carbapenem Resistance Diagnostic Kit	12 tests/kit	5-plexes
	106-0142-01	Vancomycin-Resistant Enterococcus Gene (vanA, vanB) Detection Kit	12 tests/kit	3-plexes
	106-0141-01	Staphylococcus aureus and Methicillin-Resistant Staphylococcus aureus Drug-Resistance Gene Detection Kit	12 tests/kit	3-plexes
Intestinal Infectious Diseases	106-0054-01	Coxsackievirus A6/ A10/ A16/ Enterovirus 71 RNA Diagnostic Kit	12 tests/kit	4-plexes
	106-0053-01	Coxsackievirus A16/ Enterovirus 71/ Enterovirus RNA Diagnostic Kit	12 tests/kit	3-plexes
STIs (Sexually Transmitted Infections)	106-0137-01	Chlamydia trachomatis/ Neisseria gonorrhoeae/ Ureaplasma urealyticum Diagnostic Kit	12 tests/kit	3-plexes
	106-0064-01	Nucleic acid detection kit for 16&18 human papillomavirus(HPV) genotypes	12 tests/kit	2-plexes
	106-0066-01	Diagnostic Kit for detecting of 16 high-risk and genotyping of 16&18 human papillomavirus (HPV)	12 tests/kit	5-plexes
Arbovirus Infectious Diseases	106-0044-01	Dengue virus/ Chikungunya virus/ Zika virus Diagnostic Kit	12 tests/kit	3-plexes
	106-0043-01	Dengue Serotyping Diagnostic Kit	12 tests/kit	4-plexes



Igenesis(Shanghai)Co.,Ltd.
 4/F,No.1,Lane 500,Furonghua Road,Pudong,Shanghai,P.R.China
 Tel:+86-21-38016598
<http://www.igenesisbioglobal.com>
 Mailbox:info@igenesisbio.com