



# EZER™ Adenovirus Antigen Rapid Test



## INTENDED USE

The EZER™ ADV Antigen Rapid Test is an *in vitro* diagnostic test for the qualitative detection of Adenovirus (ADV) nucleoprotein antigens in nasopharyngeal (NP) swab, throat swab and nasal aspirate samples, using the rapid immunochromatographic method.



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.



+86-571 878 18163



[www.genesis-ivd.com](http://www.genesis-ivd.com)



No.139, St.10th (East), Hangzhou Economic & Technological Development Zone, China

## PERFORMANCE CHARACTERISTICS

| Adenovirus       | Relative Sensitivity | Relative Specificity |
|------------------|----------------------|----------------------|
| Nasal Aspiration | 94.2%                | 98.7%                |
| Throat Swabbing  | 98.6%                | 99.2%                |
| NP swab          | 100%                 | 100%                 |

## AVAILABILITY

| Product         | EZER™ Adenovirus Antigen Rapid Test |
|-----------------|-------------------------------------|
| Cat. No.        | P211005                             |
| Test Device     | 20                                  |
| Sterilized Swab | 20                                  |
| Extraction Tube | 20                                  |
| Nozzle          | 20                                  |
| Package Insert  | 1                                   |


## TEST PROCEDURE

### Sample Preparation



Nasal Aspiration

OR




Throat Swabbing

OR



Nasopharyngeal swab

### Sample Preparation



## RESULTS



Positive



Negative



Invalid



Invalid