






EZER™

Strep.A Antigen Rapid Test



INTENDED USE

The EZER™ Strep A Antigen Rapid Test is a lateral flow, one-step immunoassay for the rapid, qualitative detection of *Group A streptococcus* antigen from throat swabs.



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.



+86-571 878 18163



www.genesis-ivd.com



310018



No.139, St.10th (East), Hangzhou Economic & Technological Development Zone.


PERFORMANCE CHARACTERISTICS

Strep A		Another Product		
		+	-	Total
EZER™ Strep A Antigen Rapid Test	+	116	11	127
	-	8	390	399
	Total	124	401	525



Positive Percent Agreement: 93.55%

Negative Percent Agreement: 97.26%


Overall Percent Agreement: 96.38%





RESULTS



Positive



Negative



Invalid

AVAILABILITY

Product	EZER™ Strep A Antigen Rapid Test
Cat. No.	P21116
Test Device	20
Sterilized Swab	20
Nozzle	20
Plastic Tube	20
Extraction Buffer(A/B)	2
Tube Stand	1
Package Insert	1