




*EZER*TM Dengue NS1 Rapid Test Device



INTENDED USE

The *EZER*TM Dengue NS1 Rapid Test Device is a lateral flow chromatographic immunoassay for the qualitative detection of Dengue NS1 in human whole blood, serum or plasma.

DENGUE



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.



+86-571 878 18163



www.genesis-ivd.com



310018




No.139, St.10th (East), Hangzhou Economic & Technological Development Zone.

PERFORMANCE CHARACTERISTICS


Dengue NS1 Rapid Test	Commercial EIA		
	+	-	Total
+	113	0	113
-	3	384	387
Total	116	384	500

Sensitivity: 97.4% Specificity: 100% Overall Agreement: 99.4%

PROCEDURE



RESULTS



AVAILABILITY

Product	EZER™ Dengue NS1 Rapid Test Device
Cat. No.	P231104
Test Device	20
Plastic Dropper	20
Sample Buffer	1
Package Insert	1