



# **EZER™ Dengue IgG/IgM Rapid Test Device**



## **INTENDED USE**

The *EZER™* Dengue IgG/IgM Rapid Test Device is a lateral flowchromatographic immunoassay for the qualitative detection of Dengue IgG/IgM in humanwhole blood, serum or plasma.

## **DENGUE**



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.

+86-571 878 18163

[www.genesis-ivd.com](http://www.genesis-ivd.com)


310018

No.139, St.10th (East), Hangzhou Economic & Technological Development Zone.


# PERFORMANCE CHARACTERISTICS

| EZER™ Dengue IgG/IgM Rapid Test Device | Relative Sensitivity | Relative Specificity |
|----------------------------------------|----------------------|----------------------|
| IgG/IgM vs Elisa                       | 97.40%               | 100%                 |

# PROCEDURE



# RESULTS



# AVAILABILITY

| Product                        | EZER™ Dengue IgG/IgM Rapid Test Device |
|--------------------------------|----------------------------------------|
| Cat. No.                       | P231106                                |
| Test Device                    | 20                                     |
| IgG/IgM Plastic Dropper (10µl) | 20                                     |
| Sample Buffer                  | 1                                      |
| Package Insert                 | 1                                      |