

EZER™

Treponema Pallidum Antibody Rapid Test



+86-571 878 18163



310018



www.genesis-ivd.com




No.139, St.10th (East), Hangzhou Economic
& Technological Development Zone.



INTENDED USE

The EZER™ Treponema Pallidum Antibody Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of antibodies to *Treponema Pallidum* in human whole blood/ serum/ plasma. It provides an aid in the diagnosis of syphilis.



PERFORMANCE CHARACTERISTICS



Treponema Pallidum Antibody	PCR		
	+	-	Total
+	330	1	331
-	49	630	679
Total	379	631	1010

Sensitivity: 87.1%


Specificity: 99.8%

Overall Percent Agreement: 95.0%

PROCEDURE



RESULTS



AVAILABILITY

Product	EZER™ Treponema Pallidum Antibody Rapid Test
Cat. No.	P231111
Test Device	40
Sample Buffer	1
Plastic Dropper	40
Package Insert	1