



EZER™

Strep. pneumoniae

Antigen

Rapid Test



INTENDED USE

The EZER™ *Strep. pneumoniae* Antigen Rapid Test is used for qualitative detection of *Streptococcus pneumoniae* antigen in urine to support the diagnosis of *Streptococcus pneumoniae* infection.



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.



+86-571 878 18163



www.genesis-ivd.com



310018



No.139, St.10th (East), Hangzhou Economic
& Technological Development Zone.

PERFORMANCE CHARACTERISTICS


Strep. pneumoniae Antigen Rapid Test	Another Product		
	+	-	Total
+	18	0	18
-	0	32	32
Total	18	32	50

Positive result consistency: 100.0%

Negative result consistency: 100.0%

Overall result consistency: 100.0%


PROCEDURE



AVAILABILITY

Product	<i>EZER™ Strep. pneumoniae</i> Antigen Rapid Test
Cat. No.	P241108
Test Device	40
Transfer Pipette	40
Package Insert	1


RESULTS




Positive



Negative



Invalid



Invalid