


EZER™ MCR-1 Enzyme Rapid Test

MCR-1



- ✓ High Accuracy Compared with PCR
- ✓ Fast Results in 15 Minutes
- ✓ Easy to Use, No Equipment Required
- ✓ Easy Visual Interpretation

INTENDED USE

The EZER™ MCR-1 Enzyme Rapid Test is an immunochromatographic assay (ICA) for directly qualitative detection of MCR-1 protein in bacterial culture. It's an in vitro diagnostic test for professional use only.




PERFORMANCE CHARACTERISTICS


Result		Molecular Method		
		+	-	Total
EZERTM MCR-1 Enzyme Rapid Test	+	105	0	105
	-	0	138	138
	Total	105	138	243

Relative Sensitivity: 100.00%
 Relative Specificity: 100.00%
 Accuracy: 100.00%

PROCEDURE



RESULTS



AVAILABILITY

Product	EZERTM MCR-1 Enzyme Rapid Test
Cat. No.	P261105
Test Device	20
Extraction Tube with Buffer	20
Package Insert	1