




***EZER™* Malaria Pf Antigen Rapid Test Device**



INTENDED USE

The *EZER™* Malaria Pf Antigen Rapid Test Device is a lateral flow chromatographic immunoassay for the qualitative detection of *P. falciparum* in human whole blood.

MALARIA



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.



+86-571 878 18163



www.genesis-ivd.com



310018



No.139, St.10th (East), Hangzhou Economic & Technological Development Zone.

PERFORMANCE CHARACTERISTICS

Result		Microscopy		
		+	-	Total
EZER™ Malaria Pf Antigen Rapid Test Device	+	47	0	47
	-	2	197	199
	Total	49	197	246


Relative Sensitivity: 95.92% (95%CI*: 86.02% ~ 99.50%)

Relative Specificity: >99.99% (95%CI*: 98.14% ~ 100.00%)


Accuracy: 99.19% (95%CI*: 97.09% ~ 99.90%)

PROCEDURE

Whole Blood Specimen



RESULTS



AVAILABILITY

Product	EZER™ Malaria Pf Antigen Rapid Test Device
Cat. No.	P231107
Test Device	20
Plastic Dropper	20
Sample Buffer	1
Package Insert	1